

FKG Race[®] EVO + Options



Optional
Race[®] EVO
25/.08




DIFFICULT CASES

K File
ISO 10



Race[®]EVO
Sequence .04



800 - 1000/min

Optional
XP-endo[®] Finisher
Cleaning




60 sec, 1000/min

Optional
Race[®] EVO
40/.04



LARGE CANALS

Optional
Race[®] EVO
50/.04



REFER TO FKG USER MANUAL WWW.FKG.CH

Link for videos: www.ftcdental.ch/fkg-race-evo

Obturation .02, .04, .06



**Sterile Paper
points**



**Total Fill®
BC points**



**Total Fill® BC
sealer**

REFER TO FKG USER MANUAL WWW.FKG.CH



Link for videos: www.ftcdental.ch/fkg-race-evo